




Scan
Run Program
ABORT

Scan Parameters
ID: 0
SCLN: 1000 μm
RES: 0.256 $\mu\text{m}/\text{sample}$
FORCE: 15 mg
M RANGE: 1mm
PROFILE: Hills&Valleys
LEVEL: No
R CUR: 100.00 μm
M CUR: 1000.00 μm
D RANGE: Auto

Analytic Results



R Cursor
2720 Å at 100.00 μm
M Cursor
26880 Å at 1000.00 μm
Vert D
24160 Å
Horiz_D
900.00 μm

Navigation controls: Up, Down, Left, Right arrows, F, R, M buttons.